



AviSeq™
Genomic DX

COLUNG

REF:AVG102016 16 TESTS
REF:AVG102048 48 TESTS



CONTENTS

CONTENTS	2
1. INTENDED USE	3
2. KIT CONTENTS	3
3. STORAGE AND PRODUCT STABILITY	3
4. KIT FEATURES	4
5. REQUIRED MATERIAL NOT INCLUDED	4
6. IMPORTANT NOTES AND SAFETY INFORMATION	4
7. OPERATIVE PROTOCOL	5
7.1 DNA SAMPLES PREPARATION	5
7.2 TARGET ENRICHMENT	5
7.3 LIBRARY PURIFICATION	7
7.4 QUANTIFICATION OF PCR PRODUCTS	7
7.5 QUALITATIVE AND QUANTITATIVE ANALYSIS OF LIBRARIES (OPTIONAL)	8
INDEXING/BARCODING PROTOCOL	8
7.6 INDEXING/BARCODING PCR MIX PREPARATION	8
7.7 PCR PURIFICATION	9
7.8 QUALITATIVE AND QUANTITATIVE EVALUATION OF PCR PRODUCTS (OPTIONAL)	10
7.9 QUANTITATIVE EVALUATION OF PCR PRODUCTS	10
8. TROUBLESHOOTING	11
9. SYMBOLS	11

1. INTENDED USE

Aviseq COLUNG is a kit for the analysis of the genomic targets listed in table 1. Analysis it's performed with a molecular protocol based on Next Generation Sequencing (NGS) technologies.

Aviseq COLUNG must be used together with the following AviSeq kits (not included in the kit) containing the specific adapters and index or barcode, and it's compatible with Illumina and Ion Torrent sequencers:

Illumina: Aviseq Index Set 16 (Ref. AVG404), or Aviseq Index Set 1 (Ref. AVG411)

Ion Torrent: Aviseq Barcode Set 1-16 (cod. AVG502) or Aviseq Barcode Set 17-32 (cod. AVG503)

The kit is validated for somatic and germline analysis (SNPs, indel) of DNA extracted from tissues (fresh, frozen, FFPE, FNA, etc.).

GENE	EXONS
BRAF	11, 15
PIK3CA	10, 21
EGFR	18, 19, 20, 21
KRAS	2, 3, 4
NRAS	2, 3, 4

Table 1: Detail of genes and exons targeted

2. KIT CONTENTS

The kit size is 16 or 48 tests, which can be divided into a maximum of three analysis sessions. Dividing the kits into more than three analysis sessions may alter the overall performance of the kit and / or decrease the total number of tests that can be performed.

Tube	Description	Volume 16 tests	Volume 48 tests
1	Primers mix Pool CL1	20 µl	55 µl
2	Primers mix Pool CL2	20 µl	55 µl
Taq	Hot Start Hi-Fi Taq	12 µl	32 µl
BA	BUFFER A	320 µl	920 µl
BB	BUFFER B	320 µl	920 µl
NTP	dNTPs Mix 10mM	36 µl	92 µl

Table 2: kit contents

3. STORAGE AND PRODUCT STABILITY

All reagents provided with our kits are ready to use and should be stored at -20° C.

The kit, intact and properly stored, will maintain high quality performance capacity until the expiry date marked on each single reagent tube and on the external jar/container.

When the reagents are thawed they must be kept on ice throughout the process.

4. KIT FEATURES

- SPECIFICITY: >99%
- SENSITIVITY: >99%

5. REQUIRED MATERIAL NOT INCLUDED

5.1 Generic Material

- Computer with constantly updated and guaranteed secure internet connection
- Micropipettes calibrated and periodically verified 0.2-2 µl, 2-20µl, 20-200µl or 100-1000µl and filter tips
- Vortex
- Disposable Gloves without powder
- Thermal cycler calibrated and periodically verified
- Tubes and Caps or 96-wells plate, as needed, DNase and RNase free
- Nuclease-free water
- 1,5 ml tube magnetic separator or 96-wells plate compatible magnetic separator
- Fresh 70% ethanol
- Illumina or Ion Torrent sequencers calibrated and periodically verified
- For Illumina sequencers: Aviseq Index Set 16 (Ref. AVG404), or Aviseq Index Set 1 (Ref. AVG411)
- For Ion Torrent sequencers: Barcode set 1-16 or Barcode Set 17-32 (AviSeq, cod. AVG502 or cod. AVG503)

5.2 Specific Material

The material listed below has been used and validated by AviSeq:

- Qubit™ 2.0 Fluorometer (Invitrogen Cod. Q32866) or Qubit™ 3.0 Fluorometer (Invitrogen Cod. Q33216) or Qubit™ 4.0 Fluorometer (Invitrogen Cod. Q33226) calibrated and periodically verified
- Qubit™ assay tubes (Invitrogen Cod. Q32856)
- Qubit™ dsDNA HS Assay Kit (Invitrogen, cod. Q32851)
- Magnetic beads for the purification of genomic libraries, AMPure XP (Beckman Coulter, cod. A63880)

QUALITATIVE ANALYSIS OF DNA (Optional)

- Agilent 2100 Bioanalyzer system with DNA reagent kit calibrated and periodically verified

6. IMPORTANT NOTES AND SAFETY INFORMATION

The user is required to apply the following provisions.

If the device or the results it generates, even in part, is transferred to a third party, the user must inform the end-user about the application of the specific provisions. The manufacturer is committed to constantly checking the possibilities of implementing the procedures, providing support to users.

- The kit is for professional use, it must be used by trained professionals in molecular biology.
- Do not use if package damaged.
- Biological samples and all reagents should be used in properly equipped rooms, clean and clear of potential contaminants. We suggest cleaning working areas frequently using a solution containing sodium hypochlorite 1-5%.
- Always use safety equipment such as laboratory coat, gloves and safety goggles during all steps described in the protocol.
- Check the risks and safety procedures associated with instruments, electricity, chemicals and other resources applied to the use of the device.

- Prepare ways of detecting errors in the operation of the device, evaluating after each usage the quality of the results generated; in case of doubts or anomalies found, the supplier must be promptly contacted for support.
 - When the results produced are used in diagnostic or clinical processes, the user is required to consider the possible risks associated with diagnostic errors, to set up control mechanisms and to inform the medical personnel responsible for the diagnostic or clinical processes.
 - To avoid contamination of reagents we recommend using DNase/RNase free tubes, filter tips and to pay particular attention to keep all instruments clean and free of contaminants.
 - We suggest preparing a unidirectional workflow from the initial phase of DNA isolation following the PCR preparation phase, amplification and post-amplification phases in order to keep working areas separated for the different phases of the procedure using for each phase dedicated laboratory coats, micropipettes, tubes and filter tips.
 - Used reagents and biological samples must be wasted according to legal procedures.
- **Stopping point:** every time is present a stopping point you can proceed with the following step, or store the amplified samples at 4 °C for 24 hours or -20 °C for a longer period.

7. OPERATIVE PROTOCOL

7.1 DNA SAMPLES PREPARATION

Use any commercial kit to obtain DNA from biological tissues. Determine the starting DNA concentration by fluorometric methods for accuracy. Dilute each DNA sample in PCR grade water to a final concentration of 10-25 ng/μl.

NB: The efficiency of the amplification process strongly depends on the quality of the DNA. For not fragmented DNA (i.e. fresh tissue, cytological samples...) 10 ng are in general sufficient for a good amplification process. For fragmented DNA (i.e. FFPE samples) use 25-50 ng of DNA for each PCR reaction. For highly fragmented DNA is possible to increase the quantity of DNA used up to 100-150 ng for reaction.

7.2 TARGET ENRICHMENT

Amplification of the target regions are performed with PCR in 2 reaction tubes, each containing one different primer pools. At the end of PCR, all the PCR reactions will be mixed into a single pool.

Refer to the table below for the quantities of reagents required to prepare each reaction of target enrichment:

Multiplex POOL 1 (Mix x 1 reaction)		Multiplex POOL 2 (Mix x 1 reaction)	
Primer Mix Pool CL1	1μl	Primer Mix Pool CL2	1μl
Buffer A	4μl	Buffer A	4μl
Buffer B	4μl	Buffer B	4μl
dNTPs Mix 10 mM	0,4μl	dNTPs Mix 10 mM	0,4μl
Hot Start Hi-Fi Taq	0,2μl	Hot Start Hi-Fi Taq	0,2μl
DNA	20-100ng	DNA	20-100ng
H ₂ O	Up to 20μl	H ₂ O	Up to 20μl
TOT	20μl	TOT	20μl

Table 3: quantities of reagents for each target enrichment reaction

The Reaction could be set up in 2 STEPS:

- Step 1: Prepare AmpliMIX

In a clean and appropriated work area (e.g. laminar workflow) prepare the AmpliMIX, the amplification mix containing the PCR reagents common to all reaction tubes.

- a) Unfreeze Buffer A, Buffer B, dNTPs, Primer Mix Pool 1 and 2. Vortex and briefly centrifuge. Keep Taq Polymerase in ice (do not vortex).
- b) Prepare AmpliMIX as described in the following table:

AmpliMIX	x 1
Buffer A	8 µl
Buffer B	8 µl
dNTPs Mix 10 mM	0,8 µl
Hot Start Hi-Fi Taq	0,4 µl
H ₂ O	16,8 µl
TOTAL	34 µl

Table 4: AmpliMIX for 1 sample

- Step 2: Prepare Master Mixes

- a) Accurately mix the AmpliMIX (do not vortex), then divide the volume in 2 x 1.5 ml tubes and add a specific Primer Mix Pool to each tube, as to the table below:

MMixes	x 1	
	MMix 1	MMix 2
AmpliMIX	17 µl	17 µl
Primer Mix Pool CL1	1 µl	-
Primer Mix Pool CL2	-	1 µl
TOTAL	18 µl	18 µl

Table 5: Preparation of MasterMixes

- b) For each sample, dispense 18 µl of each MasterMix per tube, change working area then add 2 µl of DNA (20-100 ng) in each of the 2 PCR tubes. Briefly centrifuge the samples.

N.B. During all time keep the mixes and reagents on ice when not in use.

Place PCR tubes in the thermocycler and then start the program according to the following thermal profile.




Table 6: PCR thermal profile

At the end of PCR, mix the PCR products deriving from the same sample in a unique tube to obtain a library pool, and then proceed with the protocol.

➤ **STOPPING POINT**

7.3 LIBRARY PURIFICATION

In appropriate post-PCR area, proceed with the library purification:

N.B. Fresh ethanol 70 % is required for a proper wash low molecular weight molecules. Low ethanol concentrations can lead to loss of sample.

- For each sample, distribute 32 µl (0.8x) of magnetic beads (previously resuspended) in low binding 1,5 ml tubes.
- Add 40 µl of each amplified library pool; pipet up and down 5 times to thoroughly mix the beads suspension with the DNA.
- Incubate the mixture for 5 minutes at room temperature.
- Place tubes on magnetic rack and incubate for 2 minutes or until the mixture appears clear; remove supernatant without touching the pellet.
- Add 200 µl of 70% Ethanol and rotate tubes twice to wash the beads; remove supernatant without touching the pellet.
- Repeat the washing step.
- Ensure that all ethanol is removed; keeping tubes on the magnetic rack let the beads dry on air for 5/10 minutes (do not desiccate).
- Remove the tubes from the magnetic rack and elute the DNA in 45µl of RNase/DNase free distilled water; pipet up and down 5 times to thoroughly mix the beads suspension with the DNA.
- Incubate 5 minutes at room temperature.
- Place tubes on magnetic rack and incubate for 2 minutes or until the mixture appears clear; carefully collect 40 µl of the supernatant containing the desired amplicons without touching the pellet.

➤ **STOPPING POINT**

7.4 QUANTIFICATION OF PCR PRODUCTS

Analyse the obtained libraries with a fluorometer Qubit and kit Qubit dsDNA HS Assay.

7.5 QUALITATIVE AND QUANTITATIVE ANALYSIS OF LIBRARIES (Optional)

A quick estimation of DNA products size may be obtained running samples on agarose gel. However, we strongly recommend using Bioanalyzer 2100 with DNA reagent kit; please refer to procedure of instrument.

INDEXING/BARCODING PROTOCOL

7.6 INDEXING/BARCODING PCR MIX PREPARATION

7.6.1 Defrost, quickly vortex and centrifuge the reagents supplied. Keep the Taq Polymerase on ice (do not vortex).

N.B. For PCR use Buffer A, Buffer B, dNTPs and Taq Polymerase contained in HotSpot panel kit; during all time keep the mixes and reagents on ice when not in use.

For Illumina sequencers: Aviseq Index Set 16 (Ref. AVG404), or Aviseq Index Set 1 (Ref. AVG411)

For Ion Torrent sequencers: Barcode set 1-16 or Barcode Set 17-32 (AviSeq, cod. AVG502 or cod. AVG503)

7.6.2 For each sample, prepare a mix according to the following scheme.

ILLUMINA	
Master Mix	x1
DNA (from purified library)	100 ng
Buffer A	10 µl
Buffer B	10 µl
dNTPs Mix 10 mM	1 µl
Hot Start Hi-Fi Taq	0.2 µl
Index i7	1 µl
Index i5	1 µl
H ₂ O	up to 50 µl
TOT	50 µl

ION TORRENT	
Master Mix	x1
DNA (from purified library)	100 ng
Buffer A	10 µl
Buffer B	10 µl
dNTPs Mix 10 mM	1 µl
Hot Start Hi-Fi Taq	0.2 µl
Barcode Oligo	2 µl
H ₂ O	up to 50 µl
TOT	50 µl

Table 7: Master mix preparation volumes

The Reaction could be set up in 2 STEPS

- **Step 1: Prepare AmpliMIX** (example for a library concentrated 10 ng/µl):

AmpliMIX	x1
Buffer A	10 µl
Buffer B	10 µl
dNTPs Mix 10 mM	1 µl
Hot Start Hi-Fi Taq	0.2 µl
H₂O	16,8 µl
TOT	38 µl

Table 8: PCR mix preparation volumes

- Step 2: Prepare Master Mixes

For Illumina: assign to each sample two indexes, by picking one of the oligo series 700 (i7) and one of the oligo series 500 (i5). **In order to choose the right combination between index i7 and index i5, please refer to the note “Technote Nextera low plex pooling guidelines Illumina” released by Illumina.**

For Ion Torrent: assign to each sample a specific barcode and add to each sample 2 µl of the selected barcode.

ILLUMINA		ION TORRENT	
Master Mix	x1	Master Mix	x1
AmpliMIX	38 µl	AmpliMIX	38 µl
Index i5	1 µl	Barcode	2 µl
Index i7	1 µl	TOT	40 µl
TOT	40 µl		

Table 9: Master Mix preparation

7.6.3 In a post-PCR area, add 10 µl (100 ng of DNA) of the purified library (7.3) Please note that the total volume of the PCR should be 50 µl.

7.6.4 Briefly centrifuge samples.

7.6.5 Keep tubes on ice until placement into the thermocycler.

7.6.6 Place the reaction tubes/plate in the thermal cycler and run the following PCR protocol:




Table 10: Thermal-cycler profile

➤ **STOPPING POINT**

7.7 PCR PURIFICATION

In appropriate post-PCR area, proceed with the library purification:

N.B. Fresh ethanol 70 % is required for a proper wash low molecular weight molecules. Low ethanol concentrations can lead to loss of sample.

- a) For each sample, distribute 40 µl (0.8x) of magnetic beads (previously resuspended) in low binding 1,5 ml tubes.
- b) Add 50 µl of each amplified library pool; pipet up and down 5 times to thoroughly mix the beads suspension with the DNA.
- c) Incubate samples at room temperature for 5 minutes.
- d) Place tubes on magnetic rack and incubate for 2 minutes or until the mixture appears clear; remove supernatant without touching the pellet.
- e) Add 200 µl of fresh ethanol 70% and rotate tubes twice to wash the beads; remove supernatant without touching the pellet.
- f) Repeat the washing step.
- g) Ensure that all ethanol is removed; keeping tubes on the magnetic rack let the beads dry on air for 5/10 minutes (do not desiccate)
- h) Remove the tubes from the magnetic rack and elute the DNA in 42µl of RNase/DNase free distilled water; pipet up and down 5 times to thoroughly mix the beads suspension with the DNA.
- i) Incubate 5 minutes at room temperature.
- j) Place tubes on magnetic rack and incubate for 2 minutes or until the mixture appears clear; carefully collect 40 µl of the supernatant containing the desired amplicons without touching the pellet.

➤ **STOPPING POINT**

7.8 QUALITATIVE AND QUANTITATIVE EVALUATION OF PCR PRODUCTS (OPTIONAL)

A quick estimation of DNA products size may be obtained running samples on agarose gel. However, we strongly recommend using Agilent 2100 Bioanalyzer (refer to manufacturers manual for protocol).

7.9 QUANTITATIVE EVALUATION OF PCR PRODUCTS

Quantify PCR products using a fluorometer, according to the manufacturer instructions, to prepare an equimolar pool of samples for sequencing.

For Illumina: dilute each library to 4nM. Create the libraries pool adding the same volume from each library (i.e. 5µl) in a new tube.

Process the libraries pool, following the instruction released by Illumina for each specific instrument. The final concentration is 12 pM for MiSeq, 1.8 pM for MiniSeq and iSeq 250 pM it is suggested to spike a 10% of the Phix reagent in the final reaction volume.











For Ion Torrent: dilute each library to 100 pM. Create the libraries pool adding the same volume from each library (i.e. 5µl) in a new tube.

Dilute to a suggested final concentration of 12pM and process in emulsion PCR (One Touch) and enrichment, according to the Ion Torrent procedure.

8. TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE	SUGGESTION
Absence of bands on agarose gel after electrophoresis	Wrong PCR thermal profile	Verify the PCR thermal profile and calibration then repeat the PCR reaction
	Mistakes in master mix preparation	Verify PCR mix components and repeat the PCR reaction
	Degraded reagents	Verify expiry date and storing conditions of the products
	Presence of inhibitors	Verify concentration and quality of DNA extracted using a spectrophotometer. If necessary, repeat DNA extraction.
	Low amount of DNA	Verify concentration and quality of DNA extracted using a spectrophotometer. If necessary, repeat DNA extraction.
Absence of bands on agarose gel after electrophoresis	Wrong PCR machine settings	Verify the PCR thermal profile and calibration then repeat the PCR reaction
Presence of fragments with low molecular weight	Mistakes in master mix preparation	Verify PCR mix components and repeat the PCR reaction

9. SYMBOLS

	<i>According to 98/79/CE Directive</i>		<i>Catalogue number</i>
	<i>In Vitro Diagnostic Medical Device</i>		<i>Batch code</i>
	<i>Expiration date</i>		<i>Temperature limitation</i>
	<i>Consult instruction for use (IFU)</i>		<i>Sufficient for n. tests</i>
	<i>Manufacturer</i>		<i>Do not use if package damaged</i>